

Air Film Calibration Board

NT - 125


Air Film Calibration Board is the newest neuronavigation tool for the Visor system, developed by ANT.

Straightforward and time-efficient coil positioning with the calibration board specially developed for the Magstim Air Film Coil

Before performing a neuronavigation session the position of the coil tracker on the coil has to be calibrated, enabling you to precisely track the position of your coil. Air Film Calibration Board allows you to place the Magstim Air Film coil fast and easy. It features a notch of the exact same shape as the coil, therefore, the coil simply needs to be placed in the notch; assuring you precise alignment of the coil center with the calibration board.

Reflecting spheres always visible

The calibration board has two legs which allow a reclined placement on a table. This means that during calibration the reflecting markers can be easily directed towards the camera. As a result the reflecting spheres are always visible without the need to hold the board.



For more information please contact our sales department at sales@ant-neuro.com or your local ANT representative.

